

# VivaLIFT!®

POWER RECLINER COLLECTION<sup>1</sup>

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER



WITH  
COUNTLESS  
POSITION  
COMBINATIONS

## URBANA<sup>1</sup>

Innovative technology meets contemporary style for a recliner that supports you and your busy life. The **Urbana Collection** by VivaLift!® Power Recliners is the perfect addition to your townhouse or condo. After a long day, experience the comfort of power lumbar, which can be adjusted at the touch of a button. The Urbana lets you fully recline with infinite, lay-flat positioning and experience maximum relaxation.



Scan here for  
more information  
on **VivaLift!®**









**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®


# URBANA<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION<sup>1</sup>

COLLECTION	URBANA <sup>1</sup>
MODEL	PLR-965M
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift 5'10" and Above
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	20"
SEAT DEPTH	20"
SEAT WIDTH	21"
TOP OF BACK TO SEAT	25"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33"
CHAIR WEIGHT	145 lbs.
WARRANTY	7-Yr. Pro-rated

## INDEPENDENT FUNCTIONS


-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



## STANDARD FEATURES

7. The USB remote adjusts all functions independently and allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

## STANDARD FABRIC: STONEWASH



GRANITE

GUNMETAL

MOSSY

1. VivaLift!<sup>®</sup> Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

(US) 800-800-8586 / (CAN) 888-570-1113  
prideliftchairs.com / pridemobility.com

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>