# INVACARE® BAR750 BARIATRIC BED





The **Invacare® BAR750 Bariatric Bed** features a heavy-duty, full-electric bed frame capable of supporting up to 750 lb, and comes standard with a battery backup with the patient in mind, this bed expands from 39" wide to a full 48" and from 80" in length to 88". Four ultra-quiet DC actuators easily reposition the head and foot sections, as well as the bed height.

#### **Features**

Heavy-duty frame design ensures added strength for client support

Split-spring design eases delivery and setup

Two locking casters on head section and two steering lock casters on foot section

Two high-low motors, one on head section and one on foot section

Pendant is designed to allow for dynamic client positioning

Quiet, smooth operation

Integrated battery backup, bed ends and half rails standard

Choose from three mattress width options

BARMATT Mattress options meet fire safety regulations 16CFR1633 and 16CFR1632

675 lb. patient weight capacity

### **HCPCS Code: E0302**







- 1. Heavy-duty frame
- 2. Split spring design
- 3. Locking casters

|--|--|

Model #	Description	
BAR750	Bariatric Bed, with 1 pair half rails	1993
Bed Packages		
BARPKG750-1633	BAR750/BARMATT39	
BARPKG750-1-1633	BAR750/BARMATT42	
BARPKG750-2-1633	BAR750/BARMATTEXP	

## Accessories

,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,		
BARMATT39	Bariatric Mattress, 700 lb. capacity	80" L x 39" W x 6.5" D
BARMATT42	Bariatric Mattress, 700 lb. capacity	80" L x 42" W x 6.5" D
BARMATT48	Bariatric Mattress, 700 lb. capacity	80" L x 48" W x 6.5" D
BARMATTEXP	Bariatric Mattress, 700 lb. capacity	80" x 39" - 48" W x 6" D (expandable)
BAR750BOLSTER	Bed extension from 80" to 88"	8" x 39" (expandable)









Shipping



Carton



Regulation



Limited

Electrical

Model #	Dimensions
BAR750	88"-96" L 15"-29" H 51" W

Sleep	Weight
Surface	Capacit
80-88" L	Safe Working
39", 42",	750 lb.*
48" W	Max. Patient

Weight Load: 440 lb. Load: 675 lb.

Dimensions 56" L 44" H 41" W

Compliance Warranty 5 Years Frame/Structure 3 Years Mechanical/

#### **Invacare Corporation** www.invacare.com

# USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 07-010 Rev. 11/18 181135



<sup>\*</sup>Includes patient, mattress and all accessories.