

Encore

Pedal Wheelchair

PEDAL
BRODA®




TIME B
New Product
Award 2019



New frame design, including adjustable wings and tall back option



Front pivot seat tilt maintains a proper position for self-propulsion




Optional dynamic rocking feature calms and creates contentment

Shown with short back, Matrix® Elan headrest and lateral support options

Patient, Facility, & Caregiver Benefits


- Increases quality of life and engagement with surroundings
- Allows for easy self-mobility and safer patient maneuverability
- Reduces risk of injury and falls with supportive stable positioning



Comfort Tension Seating® increases support as well as sitting tolerance

The Encore of Comfort, Functionality, and Style

UPHOLSTERY OPTIONS



SPECIFICATIONS

Tilt	Up to 33° (with Dynamic Rocking)
Weight Capacity	Up to 350 lbs. (114 kg)
Seat Widths	16", 18", 20", 22", 24", 26", 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm) 16", 18", 20", 22" with Dynamic Rocking (40cm, 45cm, 50cm, 55cm)
Seat Depths	17", 18.5", 20" (43cm, 46cm, 50cm) 15.5", 17", 18.5", 20" with Dynamic Rocking (39cm, 43cm, 46cm, 50cm)
Back Lengths	22", 32", 33.5", 35" (55 cm, 81cm, 85cm, 89cm)
Seat Surface Heights	Adjustable 12.5" - 17.5" (31cm-44cm) Adjustable 14" - 19" with Dynamic Rocking (35-48cm)
Wheels	5" Front Casters: Total Lock (12.5cm) 20" or 22" Front Mag (50cm or 55cm) 5" Rear Casters (12.5cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Tall back: Neckrest, shoulder bolsters, lower side pads, seat and back pad with cover Short back: lower side pads, seat and back pad with cover

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps