

Synthesis

Positioning Wheelchair

POSITIONING
BRODA®



Tilt-in-space functionality provides additional pressure offloading



Front pivot seat tilt maintains a forward line of sight, enhancing engagement.



Adjustable upper lateral supports fit a variety of users




Optional WC-19 Transportation Package enables safe patient travel

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Promotes safe patient handling with enhanced positioning
- Proper seating alignment reduces patient fatigue and discomfort

*The Synthesis of Comfort,
Functionality, and Style.*

UPHOLSTERY OPTIONS



SPECIFICATIONS

Tilt - Infinitely Adjustable	up to 40°
Recline - Infinitely Adjustable	up to 90°
Weight Capacity	350 lbs. (158 kg)
Seat Widths	16", 18", 20", 22", 24", 26" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm)
Seat Depths	17", 18.5", 20" (43cm, 46cm, 50cm)
Floor-to-Seat Heights	19", 20.5", 22", 23.5" (48cm, 51cm, 55cm, 59cm)
Wheels	14", 20" Rear Mount Mag (35cm, 50cm, 55cm) 14" Mid Mount Mag (35cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding	Neckrest, upper laterals, side pads, sole and calf pads, seat and back pads

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps