

THE MOST COMFORTABLE RECLINER IN THE WORLD!

Wireless Charging Station*





ATI AS PLUS

When you need a moment for yourself, the Atlas Plus Collection is there for you. Customize your comfort while you watch the big game or read your favorite book. Enjoy a refreshing beverage and keep it close with the cup holder, and conveniently charge your phone with the wireless charger. Take advantage of those rare moments of solitude with the Atlas Plus Recliner.

* The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.

**This product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Call your local Pride provider for further information. Other conditions apply.



Scan here for more information on VivaLift!®





ATLAS PLUS



ATLAS PLUS ¹
PLR-2985M
Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9"
Medium
400 lbs.
19.5″
21″
22″
26″
Yes
Standard
Total Comfort
42″
16″
39″
164 lbs.
7-Yr. Pro-rated

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/-1''. Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply. **This product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Call your local Pride provider for further information. Other conditions apply.

INDEPENDENT FUNCTIONS

- 1. The power headrest lifts your head, neck and shoulders for maximum support.
- 2. The upward lift offers convenient stand assist for getting in and out of the recliner.
- 3. The power lumbar enables you to personalize your comfort.
- 4. The footrest extension allows you to relax and stretch out completely.
- 5. The backrest moves up and down independently to help you find the perfect position.

6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

STANDARD FEATURES

- 7. The USB remote adjusts all functions independently and allows you to charge your smart device.
- 8. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: BADLANDS



MUSHROOM

STEEL

(US) 800-800-8586 / (CAN) 888-570-1113 prideliftchairs.com / pridemobility.com

